Little Doctor

LD23, LD23A, LD23L

EN INSTRUCTION MANUAL FOR DIGITAL BLOOD PRESSURE MONITOR LD23, LD23A, LD23L
PL INSTRUKCJA OBSŁUGI CIŚNIENIOMIERZA CYFROWEGO LD23, LD23A, LD23L
HU HASZNÁLATI ÚTMUTATÓ. LD23, LD23A, LD23L készülék a vérnyomás és pulzusszám digitális mérése
RO MANUAL DE UTILIZARE A TENSIOMETRULUI DIGITAL LD23, LD23A, LD23L
BG РЪКОВОДСТВО ЗА ЕКСПЛУАТАЦИЯ НА ЦИФРОВ ТОНОМЕТЪР LD23, LD23A, LD23L

fig.1 rys. 1 1. ábra
PARTS AND COMPONENTS
PODSTAWOWE CZĘŚCI I KOMPONENTY
ALKATRÉSZEK MEGNEVEZÉSE

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
fig.5 rys. 5 5. ábra
MEASUREMENT PROCEDURE
PROCEDURA POMIARU
A MÉRÉS MÓDSZERE

EN
PARTS AND COMPONENTS (fig.1)
● ELECTRONIC UNIT.
● LCD.
● ARM CUFF JACK.
● AIR PLUG.
● AIR TUBE.
● ARM CUFF.
● POWER SOURCE JACK.
● ELECTRICAL POWER SUPPLY LD-N057 (IT IS ATTACHED TO MODIFICATION LD23A AND LD23L).
● (C) (POWER ON/OFF).
● POWER ELEMENTS.
● INSTRUCTION MANUAL.
● WARRANTY CARD.
● STORAGE CASE.
GENERAL INFORMATION
This instruction Manual is designed to assist the user with safe and effective operation of the automatic digital Device for measurement of blood pressure and heartbeat rate LD (modification LD23, LD23A, LD23L) (hereinafter – the "Device"). Use the Device according to the rules described in this Manual. Operate the Device only as intended. Do not use the Device for any other purposes. Read and understand the whole Instruction Manual, in particular "Recommendations on Correct Measurement".
INDICATIONS FOR USE
Use this Device to measure your systolic and diastolic blood pressure and heartbeat rate in patients aged from 15. This Device is recommended for use by persons with unstable blood pressure or known arterial hypertension at home as an addition to medical surveillance.
OPERATION PRINCIPLE
This Device uses the oscillometric method of blood pressure with Fuzzy Algorithm and pulse rate measurement. Wrap the cuff around your upper arm and it starts to be inflated automatically. The sensitive element of the Device feels the weak pressure oscillations in the cuff generated by widening and contraction of the brachial artery in response to every heartbeat. Pumping in is ceased when cuff is adequately pumped in to determine diastolic and systolic pressure (the amplitude of the pressure waves is measured, converted into millimeters Hg and shown on the display as figures) after which air is released from cuff. Remember that the Device will not maintain the mentioned accuracy of a measurement if it is used or stored at a temperature or humidity other than those specified in Technical Specifications of this Manual. We are warning about possibility of mistakes in blood pressure measurement with this Device in persons with pronounced cardiac arrhythmia. Consult the doctor concerning blood pressure measurement of your child.

Measurement method	oscillometric with Fuzzy Algorithm
Display	LCD
Pressure indication range in an arm cuff, mmHg	from 0 to 300
Measurement range: pressure in an arm cuff, mmHg	from 40 to 260
pulse rate, 1/min	from 40 to 160
Range of admissible absolute error at measurement of air pressure in an arm cuff, mmHg	±3
Range of admissible relative error at pulse rate measurement, %	±5
Inflation	automatic (air pump)
Deflation	automatic
Type of power supply:	4 "AA" size batteries (LR6) or adapter not less than 600 mA
Operation conditions	
Temperature, °C	from 10 to 40
Relative humidity, % Rh	85 and lower
Storage and transportation conditions	
Temperature, °C	from -20 to 50
Relative humidity, % Rh	85 and lower
Date and time	yes
Max power intake, W	3,6
Output voltage	6V ± 5%
Max load current	not less than 600 mA
Input voltage	~200-240 V, 50/60 Hz
Dimensions	64 x 70 x 43 mm
Weight	not more 0,3 kg
Plug:	
Polarity of terminals	↔ internal
Internal diameter, mm	2,1 ± 0,1
External diameter, mm	5,5 ± 0,1
Length of plug contact, mm	10 ± 0,5
Dimensions	
Size (elektronic block), mm	129 x 68 x 54
Weight (without package, case, batteries and adapter), g	328 (LD23, LD23A) / 358 (LD23L)
Year of manufacture	You the manufacturer is given in the bottom of the unit body in a serial number after symbols "AA"
Opisanie simboli	± Type BF □ Important: Read the instructions Type II

Revision date of the present Manual is indicated on the last page as EXXX/YYMM/NN, where YY is the year, MM is the month and NN is the number of revision. Technical characteristics may be changed without preliminary notification to improve the operation and quality of the product.
This Device manufacturing is certified according to international standard ISO 13485. Units LD23, LD23A, LD23L comply with the requirements of European Directive MDD 93/42/EEC, international standards, EN980.
Power source LD-N057 complies with international standard EN 55022 Class A, protection level against electric shock: Class II, Type BF.
RECOMMENDATIONS ON CORRECT MEASUREMENTS
1. For correct measurement you should know that the BLOOD PRESSURE IS SUBJECT TO SHARP VARIATIONS EVEN WITHIN THE SHORT TIME INTERVALS. The blood pressure depends on many factors. It is usually lower in summer and higher in winter. The blood pressure varies together with the atmospheric pressure, depends on physical loads, emotional excitement, stresses and dietary regime. Drugs, drinking alcohol and smoking produce significant effect. Even the very procedure of blood pressure measurement in a polyclinic does not influence the blood pressure high in many people, thus, the blood pressure measured at home often differs from the values received in a polyclinic. As the blood pressure tends to rise at low temperatures, make measurements at an indoor temperature (approximately 20° C). If this Device stayed under a low temperature, keep it for at least 1 hour at an indoor temperature before use, otherwise the measurement result may be incorrect. During a day the difference in readings for healthy people may be 30-50 mmHg of systolic pressure and 10 mmHg of diastolic pressure. The dependence of the blood pressure on various factors is individual for each person. Accordingly, it is recommended to keep a special book with blood pressure records. ONLY A CERTIFIED DOCTOR USING YOUR RECORDS IS CAPABLE TO ANALYZE THE TENDENCY OF YOUR BLOOD PRESSURE VARIATIONS.
2. At cardiovascular and some other diseases requiring blood pressure monitoring make measurements in the hours fixed by your attending doctor. REMEMBER THAT THE DIAGNOSTIC AND ANY TREATMENT OF HYPERTENSION MAY BE CONDUCTED ONLY BY A CERTIFIED DOCTOR ON THE BASIS OF BLOOD PRESSURE VALUES OBTAINED BY THIS DOCTOR. TAKING OF DRUGS AND THEIR DOSES SHOULD BE PRESCRIBED ONLY BY YOUR ATTENDING DOCTOR.
3. At such disorders as deep vascular sclerosis, weak pulse wave and also in patients with the prominent distortions of cardiac rhythm it may be difficult to measure the blood pressure accurately. IN SUCH CASES CONSULT A CERTIFIED DOCTOR ABOUT APPLICATION OF THE ELECTRONIC DEVICE.
4. KEEP CUFF DURING A MEASUREMENT TO OBTAIN THE ACCURATE VALUES OF YOUR BLOOD PRESSURE WITH THE ELECTRONIC DEVICE. Measure your blood pressure in the calm and comfortable conditions at the indoor temperature. No eating an hour before measurement; no smoking, taking tonic agents, alcohol 1.5-2 hours before measurement.
5. The accuracy of blood pressure measurement depends on whether the cuff matches the size of your arm. The CUFF SHOULD NOT BE TOO SMALL OR TOO LARGE.
6. Wait 3 minutes between measurements for the blood to restore its circulation. However, the persons with prominent atherosclerosis due to considerable loss of vascular elasticity may need to increase the wait time between measurements (10-15 minutes). This also refers to the patients suffering from long from diabetes. For more accurate determination of blood pressure it is recommended to make a series of 3 consecutive measurements and to use the average value.

POWER SUPPLY OF THE DEVICE
BATTERY INSTALLATION (fig.2)
1. Open the cover of the battery compartment and install 4 "AA" size batteries according to polarity marked inside the compartment. Do not use much force to remove the cover of the battery compartment.
2. Close the battery cover.
The batteries supplied with the Device are intended for check of the Device performance at sale and their service life may be shorter than of the recommended batteries.

• Replace all batteries when the Low Battery Indicator () appears on the screen or when there is no any indication on the screen. The Low Battery Indicator does not waste the discharge level.
• Replace all four batteries at the same time. Do not use the waste batteries.
• If the Device is unused for a long time, remove all batteries.
• Do not leave the waste batteries in the Device.

USE OF THE DEVICE WITH THE POWER SOURCE
The jack for connection of the power source is on the right side of the Unit.
Use the power source only with the following technical characteristics:
The output voltage : 6V ± 5%
Maximum output current : not less than 600 mA
Plug:
Terminal polarities : "minus" – internal
External diameter : 5,5 ± 0,1 mm
Internal diameter : 2,1 ± 0,1 mm
Length : 10 ± 0,3 mm
The manufacturer recommends application of the stabilized power source LD-N057 (it is attached to modification LD23A and LD23L).
CORRECT POSITION DURING MEASUREMENT (fig.3)
1. Sit at a table so that during blood pressure measurement your hand rests on its surface. Be sure that the cuff is placed approximately at the level of your heart and that your arm lies freely on the table and does not move.

3. You can measure the pressure lying on the back. Look at the ceiling, keep your head and do not move during measurement. Be sure that the cuff is placed approximately at the level of your heart.
CUFF PREPARATION (fig.4)
1. Insert the cuff end for about 5 cm into a metal ring as shown in the figure.
2. Apply the cuff to your left upper arm so that the air tube is directed to your palm. If the measurement on your left arm is difficult, you may use your right arm. In this case remember that the readings may differ by 5-10 mmHg and even more.
3. Wrap the cuff around your upper arm so that the bottom of the cuff is approximately 2-3 cm above your elbow. The sign "ARTERY" should be over the arm artery.
4. Fix the cuff so that it fits tightly to the arm, but see that it is not overnight. Too tight or too free placement of the cuff may give inaccurate readings.
5. On the fixed cuff the sign "INDEX" should point to the area «NORMAL». It means that the cuff is chosen correctly and fits the size of your upper arm. If the sign points to the area marked «**LOW**» the cuff is too small and the readings will be higher. If the sign points to the area marked «**HIGH**» the cuff is too large and the readings will be lower.
6. If the arm has a conic form, the cuff should be put on with a spiral movement as shown in the figure.
7. If the rolled-up sleeve squeezes the arm interfering with free blood flow the Device may give inaccurate figures not corresponding to your actual blood pressure.

MEASUREMENT PROCEDURE (fig.5)
1. Insert the Air Tube Tube into the Cuff Jack. Make 3-5 deep inhaleds and exhales before taking a measurement and relax. Do not move, do not speak and do not touchen your arm.
2. Press ().
3. All symbols will appear on the display screen for a short time, two short sound signals will be given and the Device will inflate automatically the air into the cuff. At first the inflation will stop at the level of 190 mmHg.
4. After reaching the level of 190 mmHg the cuff will gradually deflate. The figures on the screen will count back. The pulse symbol ♥ will start flickering.
As the blood pressure and pulse are measured during air deflation from the arm cuff keep quiet and do not move your arm and do not touchen your arm muscles.
5. When the measurement is complete the sound signal is given, the arm cuff completely deflates and you measurement results flash on the screen and Indicator «**PA**» will blink, reminding that to retain results, it is necessary to choose memory 1 or 2, having depressed M1 or M2, respectively. The result may be before the beginning of the next measurement or before turning the device OFF.
If irregular pulse rhythm is detected during measurement, symbol of arrhythmia ♥♥♥, will appear upon measurement end (for LD23, LD23A). Appearance of arrhythmia indicator may also be caused by body movement during measurement. During periodical appearance of this indication apply to you attending doctor (for LD23, LD23A).
Apart from numerical value of pressure, result is also displayed on scale WHO (fig. 14). Scale WHO – three-color scale of classification of received value of arterial pressure, according to recommendation of World Health Organization. The scale is available from the left.
6. Press () – to switch off. For taking a new measurement repeat all steps described in this paragraph.
TO OBTAIN THE INTERNAL BETWEEN MEASUREMENTS TO RESTORE THE BLOOD CIRCULATION, WAIT FOR AT LEAST 3 MINUTES BEFORE MAKING A NEW MEASUREMENT.
If the Unit is ON and is unused for 3 minutes it will be switched off automatically.

AUTOMATIC RE-INFLATION
When during the first blood pressure measurement the cuff inflation to a level of 190 mmHg is not sufficient or you move your arm the Device stops measurement and re-inflates the cuff to the higher level. The Device has 4 fixed levels of the arm cuff inflation: 190, 230, 270 and 300 mmHg.
The automatic re-inflation is repeated until the measurement is completed successfully. This is not a defect.

FORCED DEFLATION FROM A CUFF
For rapid air release from of the arm cuff during arm cuff inflation or during a measurement (slow deflation) press the () button. The arm cuff quickly completely deflates and the Device is independently switched off.

Indication	Likely cause	Methods of correction
	The arm cuff is applied incorrectly or the air tube plus is inserted not tightly enough. Measurements cannot be made due to hand movement or talking during measurements. At such disorders as deep vascular sclerosis, weak pulse wave and also in patients with the prominent distortions of cardiac rhythm it may be difficult to measure the blood pressure accurately.	Be sure that the arm cuff is applied correctly and the tube plus is inserted tightly. Repeat the whole measurement procedure. Repeat the measurement following strictly the recommendations of this Manual. In such cases consult a certified doctor about application of the electronic unit.
	Batteries are discharged.	Replace the batteries for new ones.

CARE, STORAGE, REPAIR AND DISPOSAL
1. Keep this Device from exposure to higher humidity, direct sunlight, shocks, vibration. THIS DEVICE IS NOT WATERPROOF!
2. Do not keep and use this Device near heating installations and open fire.
3. Remove the batteries if the Device will be unused for a long time. Battery leakage may damage the Device. KEEP BATTERIES OUT OF REACH OF CHILDREN!
4. Keep the Device clean and protect it from dust. Use the dry soft cloth to clean the Device.
5. Keep the Device and its components away from water, solvents, spirit and petrol.
6. Protect the arm cuff from contacting on sharp things; do not stretch or fold tightly the arm cuff.
7. Do not subject the Device to strong shocks, such as dropping on the floor.
8. This Device does not contain special controls to adjust the measurement accuracy. It is prohibited to open individually the electronic block. Repair the Device only in authorized organizations.
9. On expiration of the warranted service life apply from time to time to authorized repair organizations to check the technical condition of the Device.
10. Dispose of the Device and its components according to the application local regulations. No special requirements to disposal of this Device are defined by the manufacturer.
11. The arm cuff may withstand multiple sanitary treatments. The internal tissue surface of the arm cuff (contacting on arm) may be cleaned with cotton ball soaked in 3%-solution of hydrogen peroxide. At long use the partial color fading of the tissue coating of the arm cuff is possible. Washing and ironing of the arm cuff are not allowed.

PROBLEM	LIKELY CAUSE	METHOD OF CORRECTION
After pressing the () button no indication on the display.	Discharge of batteries. Polarity of batteries is not observed. Power terminals are contaminated Power source not plugged in an electrical outlet.	Replace all batteries for new ones. Install batteries correctly. Clean the terminals with dry cloth. Plug the power source into an electrical outlet.
Inflation is stopped and resumed.	Automatic re-inflation to obtain correct measurements. Perhaps you talk or move your arm during the measurement?	See MEASUREMENT PROCEDURE Calm down and repeat the measurement.
Every time the blood pressure is different. Measurements are too low/high	Check that the arm cuff is level with your heart? Check that the arm cuff is applied correctly? Perhaps your arm muscles are tough? Perhaps you talk or move your arm during the measurement?	Take the correct position for measurement. Take the correct position for measurement. Calm down, apply the arm cuff correctly. Keep silence and quiet during measurement.
Measurements of the pulse rate are too high/low...	Perhaps you talk or move your arm during the measurement? Perhaps you make measurement directly after physical load?	Keep silence and quiet during measurement. Repeat the measurement at least in 5 min.
Impossibility to make a large number of measurements.	Application of poor batteries.	Use only alkali batteries of well-known manufacturers.
Spontaneous failure of power supply.	Actuation of automatic de-energizing system.	This is not a defect. The Unit is disconnected automatically in 3 minutes after the last operation of the Unit.

If regardless of the above recommendation you are unable to get correct measurement results, stop the use of the Device and apply to a maintenance organization (addresses and telephones of authorized organizations may be found in the warranty card). Do not try to adjust the internal mechanism by yourself.

WARRANTY
The warranty liabilities are contained in the warranty card given at the sale of this Device to a purchaser.
The address of organizations for warranty maintenance are given in the warranty card.
Manufactured under control and for Little Doctor International (S) Pte. Ltd. (35 Selegie Road #09-02 Parklane Shopping Mall, Singapore 188307, postal address: Yishun Central P.O. Box 9293 Singapore 917699).
Manufacturer: Little Doctor Electronic (Nantong) Co. Ltd. (No. 8, Tongjing Road Economic & Technical Development Area, 226010 Nantong, Jiangsu, People's Republic of China).
Distributor in European Union: Little Doctor Europe Sp. z o.o. (57C Zawila Street, 30-390 Kraków, Poland, phone: +48 12 2684746, fax: +48 12 268 47 53, e-mail: biuro@littledoctor.pl)
For more information please visit www.littledoctor.sg

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