







EN INSTRUCTION MANUAL FOR DIGITAL BLOOD PRESSURE MONITOR LD30

PL INSTRUKCJA OBSŁUGI CIŚNIENIOMIERZA CYFROWEGO LD30

HU HASZNÁLATI ÚTMUTATÓ. LD30 KÉSZÜLÉK A VÉNYOMÁS ÉS PULSUZZÁM DIGITÁLIS MÉRÉSÉRE

RO MANUAL DE UTILIZARE A TENSIOMETRULUI DIGITAL LD30

BG РЪКОВОДСТВО ЗА ЕКСПЛОАТАЦИЯ НА ЦИФРОВИ ТОНOMETРЪЛ LD30

fig.1 rys. 1 1. ábra imaginea 1 рис. 1  
PARTS AND COMPONENTS DENUMIREA PARTILOR SI COMPONENTELOR НАЗВАНИЯ НА ЧАСТИ И КОМПОНЕНТЫ  
ALKATRÉSZK MEGNEVEZÉSEfig.2 rys. 2 2. ábra imaginea 2 рис. 2  
BATTERY INSTALLATION INSTALAREA BATERIILOR ПОСТАВЯНИЕ НА ЕЛЕМЕНТИТЕ НА ЗАХРАНВАНЕfig.3 rys. 3 3. ábra imaginea 3 рис. 3  
CORRECT POSITION DURING MEASUREMENT ПРЫЗИЕ ПОЗИЦИИ УМОЛЖИВАЯЩЕЙ ПОМИАР А В ГЕЛЕСЫ TESTHYZETZET A MÉRÉSNEL ПОЗИТА CORRECTА IN TIMPLU PROCESULUI DE MÁSURARE ПРАВИЛНА ПОЗА ПРИ ИЗМЕРВАНИИfig.4 rys. 4 4. ábra imaginea 4 рис. 4  
CUFF PREPARATION ПРЕПАРАЦИЯ МАНКЕТУ А МАНДЗСЕТТА ЕЛОКЕШИТЕfig.5 rys. 5 5. ábra imaginea 5 рис. 5  
MEASUREMENT PROCEDURE ПРОЦЕДУРА ПОМИАРЫ НАЧИН НА ИЗМЕРВАНЕfig.6 rys. 6 6. ábra imaginea 6 рис. 6  
MEMORY FUNCTION ФУНКЦИЯ ПАМИЕЦИ А МЕРЭС МÓДСЗЕРЕ

## EN

PARTS AND COMPONENTS (fig.1)

- ① ELECTRONIC UNIT.
- ② ARM CUFF JACK.
- ③ LCD.
- ④ AIR TUBE WITH AIR PLUG.
- ⑤ Arm cuff LD.
- ⑥ Button 1 (Power ON/OFF).
- ⑦ Button M1 (memory 1).
- ⑧ Bag.

## GENERAL INFORMATION

This Instruction Manual is designed to assist the user with safe and effective operation of the automatic digital Device for measurement of blood pressure and heartbeats rate LD (modification LD30) (hereinafter - "the Device"). Use this Device according to the rules described in this Manual. Operate the Device only as intended. Do not use the Device for any other purposes. Read and understand the whole Instruction Manual, in particular "Recommendations on Correct Measurement".

## INDICATIONS FOR USE

Use this Device to measure your systolic and diastolic blood pressure and heartbeats rate in patients aged from 15. This Device is recommended for use by persons with unstable blood pressure or known arterial hypertension at home as an addition to medical surveillance. The cuff is designed for the upper arm with the circumference approximately from 25 to 36 cm.

## OPERATION PRINCIPLE

This Device uses the oscillometric method of blood pressure and pulse rate measurement. Wrap the cuff around your upper arm and it starts to be inflated automatically. The sensitive element of the Device feels the weak pressure oscillations in the cuff generated by widening and contraction of the brachial artery in response to every heartbeat. Pumping is ceased when cuff is adequately pumped up to determine diastolic and systolic pressure (the amplitude of the pressure waves is measured, converted into millimeters Hg and shown on the display as figure) after which air is released from cuff. The Device has 2 memories, by 60, for storage of measurement results. Remember that the Device will not maintain its mentioned accuracy if it is used stored at a temperature or humidity other than those specified in Technical Specifications of this Manual. We are warning about possibility of mistakes in blood pressure measurement with this Device in persons with pronounced cardiac arrhythmia. Consult the doctor concerning blood pressure measurement of this Device.

The pulse rhythm is detected during measurement, symbol of arrhythmia (♥), will appear upon measurement end. Appearance of arrhythmia indication may also be caused by body movement during measurement. During periodical appearance of arrhythmia pulse rhythm is detected during measurement.

Apart from standard value of pressure result is also displayed scale WHO (Fig. 14). Scale WHO - three-color scale of classification of received value of arterial pressure, according to recommendation of World Health Organization. The scale is available from the left. Press button ☰ to shift to mode of date and time setting. TO OBTAIN THE ACCURATE RESULT MAKE INTERVAL BETWEEN MEASUREMENTS AS SHORT AS POSSIBLE. WAIT FOR AT LEAST 3 MINUTES BEFORE MAKING A NEW MEASUREMENT.

THE DATA WILL NOT BE KEPT IN THE MEMORY EVEN IF THE DEVICE IS STORED WITHOUT BATTERIES. TO DELETE ALL VALUES STORED IN THE MEMORY YOU SHOULD MAKE ACTIONS DESCRIBED IN "MEMORY FUNCTION".

If no actions are made in mode of date and time setting for more than 3 minute, the Device is independently changed over to mode of date and time indication.

## AUTOMATIC RE-INFLATION

When during the first blood pressure measurement the cuff inflation to a level of 190 mmHg is not sufficient or you move your arm the Device stops measurement and re-inflates the cuff to the higher level. The Device has 4 fixed levels of the arm cuff inflation: 190, 230, 270 and 300 mmHg.

The automatic re-inflation is repeated until the measurement is completed successfully. This is not a defect.

## FORCED DEFLECTION FROM A CUFF

For rapid air release from the arm cuff during cuff inflation or during a measurement (slow deflation) press the ☰ button. The device will quickly release all air from cuff and will change over to mode of date and time indication.

## MEMORY FUNCTION (fig.6)

1. Result of each measurement (pressure, pulse, time and date of measurement) may be kept in the device memory. For this purpose, after measurement, within not more than 3 minutes, memory M1 or M2 shall be selected for memory storage. IF THE NOTICE ON ERROR APPEARS THE MEASUREMENT RESULT WILL NOT BE STORED.

2. Up to 60 measurement results and mean value of 3 measurements may be kept in each device memory. When the number of measurements exceeds 60, the oldest record is deleted to save the most recent values.

3. Press the M1 (or M2) button to see the figures stored in the memory. At the first depression of button M1 (or M2) mean value of 3 last measurements, kept in memory M1 (or M2) and number of memory cell will appear on the screen, and in next 3 seconds its content is displayed (Fig. 6). When content of memory cell is displayed, date and time of measurement are displayed alternately in the display lower line. Each depression of button M1 (or M2) causes shifting to the next memory cell.

4. To delete all values stored in the memory, press the M1 (or M2) button and hold it down for more than 5 seconds. Symbols "C" will appear on the screen and all values will be deleted from the memory.

## INFORMATION ABOUT ERRORS

**Err** The arm cuff is applied incorrectly or the air tube plus is inserted not tightly enough.

Measurements cannot be made due to hand movement or talking during measurements.

Repeat the measurement following strictly the recommendations of this manual.

**Err** Batteries are discharged.

Replace the batteries for new ones.

**Care, Storage, Repair and Disposal**

1. Keep this Device from exposure to high humidity, direct sunlight, shocks, vibration. THIS DEVICE IS NOT WATERTIGHT!

2. Do not keep and use this Device near heating installations and open fire.

3. If the Device was stored at a temperature below the freezing point, keep it at least for 1 hour in some warm place before use.

4. Remove the battery if the Device will be unused for a long time. Battery leaking may damage the Device. KEEP BATTERIES OUT OF REACH OF CHILDREN!

5. Keep the Device clean and protect it from dust. Use the dry soft cloth to clean the Device.

6. Keep the Device and its components away from water, solvent, spirit and petrol.

7. Protect the arm cuff from contacting on sharp things; do not stretch or fold tightly the arm cuff.

8. Do not subject the Device to strong shocks, such as dropping on the floor.

9. This Device does not contain special controls to adjust the measurement accuracy. It is prohibited to open individually the electronic block. Repair the Device only in authorized organizations.

10. On expiration of the warranty service life apply from time to time to authorized repair organizations to check the technical condition of the Device.

11. Dispose of the Device and its components according to the application local regulations. No special requirements to disposal of this Device are defined by the manufacturer.

12. The arm cuff may withstand multiple sanitary treatments. The internal tissue surface of the arm cuff (contacting a skin) may be cleaned with cotton ball soaked in 3% solution of hydrogen peroxide. At long use the partial color fading of the cuff is possible. The dependence of the blood pressure on various factors is individual for each person. Accordingly, it is recommended to keep a special book with pressure values. ONLY A CERTIFIED DOCTOR USING YOUR RECORDS IS CAPABLE TO ANALYZE THE TENDENCY OF YOUR BLOOD PRESSURE VALUES.

13. Do not use vascular and somatic diseases monitoring blood pressure monitoring mode measurements in the hours fixed by your attending doctor. REMEMBER THAT THE DIAGNOSTIC AND ANY TREATMENT OF HYPERTENSION MAY BE CONDUCTED ONLY BY A CERTIFIED DOCTOR ON THE BASIS OF BLOOD PRESSURE VALUES OBTAINED BY THIS DOCTOR. TAKING OF DRUGS AND THEIR DOSES SHOULD BE PREScribed ONLY BY YOUR ATTENDING DOCTOR.

14. Do not use the Device as a drug for the prevention of stroke pulse wave and also in patients with the prominent distortions of the vessels coating of the cuff is possible. Washing and ironing of the arm cuff are not allowed.

15. Do not use the Device for the blood to restore its circulation. However, the persons with prominent arteriosclerosis due to considerable loss of vascular elasticity may need to increase the wait time between measurements (10-15 minutes). This also refers to the patients suffering from long diabetes. For more accurate determination of blood pressure it is recommended to make a series of 3 consecutive measurements and to use the average value.

**POWER SUPPLY OF THE DEVICE**

BATTERY INSTALLATION (fig.2)

1. Open the cover of the battery compartment and install 4 "AA" size batteries according to polarity marked inside the compartment. Do not use much force to remove the cover of the battery compartment.

2. Close the battery cover.

The batteries supplied with the Device are intended for check of the Device performance at sale and their service life may be shorter than the recommended batteries.

• Replace all batteries when the Low Battery Indicator **■** appears on the screen or when there is no indication on the screen. The Low Battery Indicator does not show the discharge level.

• If the Device is unused for a long time, remove all batteries.

• Do not leave the waste batteries in the Device.

USE OF THE DEVICE WITH THE POWER SOURCE

The manufacturer recommends application of the stabilized power source LD-N05T (it is attached to modification LD30). The jack for connection of the power source is on the right side of the Device.

The addresses of organizations providing warranty service can be found in section "INFORMATION ON THE MANUFACTURER AND DISTRIBUTORS" or on the website www.littledoctor-sg.

INFORMATION ON THE MANUFACTURER AND DISTRIBUTORS

Manufactured under control and for Little Doctor International (S) Pte. Ltd. (7500A, Beach Road, 11-13 The Plaza 199591, Singapore. Postal address: Yishun Central P.O. Box 9293 Singapore 917699).

Manufacturer: Little Doctor Electronic (Nantong) Co., Ltd. (No. 8, Tongxing Road Economic & Technical Development Area, 226010 Nantong, Jiangsu, People's Republic of China).

Distributor in European Union: Little Doctor Europe Sp. z o.o. (576 Zawila Street, 30-390 Krakow, Poland, phone: +48 12 2684746, 12 2684747, fax: +48 12 268 47 53, e-mail: biuro@littledoctor.pl)

For more information please visit www.littledoctor.sg

CORRECT POSITION DURING MEASUREMENT (fig.3)

1. Sit at a table so that during blood pressure measurement your hand rests on its surface. Be sure that the cuff is placed approximately at the level of your heart and that your arm lies freely on the table and does not move.

2. You can measure the blood pressure lying on the back. Lie on the ceiling, keep quiet and do not move during measurement. Be sure that the cuff is placed approximately at the level of your heart.

CUFF PREPARATION (fig.4)

1. Insert the cuff end for 5 cm into a metal ring as shown in the figure.

2. Apply the cuff to your left upper arm so that the air tube is directed to your palm. If the measurement on your left arm is difficult, you may use your right arm. In this case remember that the readings may differ by 10-15 mmHg and even more.

3. Wrap the cuff around your upper arm so that the bottom of the cuff is approximately 2-3 cm above your elbow. The sign "ARTERY" should be over the arm artery.

4. Fix the cuff so that it fits tightly to the arm, but see that it is not overtight. Too tight or too free placement of the cuff may give inaccurate readings.

5. On the fixed cuff the sign "INDEX" should point to the area "NORMAL" (25-36 cm). It means that the cuff is chosen correctly and fits the size of your upper arm. If the sign points to the area marked **■**, the cuff is too large and the readings will be lower.

6. If the arm has a conic form, the cuff should be put with a slight movement as shown in the figure.

7. If the rolled-up sleeve interferes with free blood flow the Device may give inaccurate figures not corresponding to your actual blood pressure.

SETTING OF DATE AND TIME

To shift mode of date and time setting, it is necessary, holding button M1, depress button. Selected parameter will blink. Change of selected parameter is taking place when depressing button M2 towards decreasing - when depressing button M1. To shift to setting of the next parameter, year/month/year, it is necessary to depress. If action done in mode of date and time setting, indication of date and time will stop. Measurement of pressure and pulse is possible without date and time setting. Indication of date and time is switched off when device is used as clock, one set of power cells will be enough approximately for 1 year.

MEASUREMENT PROCEDURE (fig.5)

1. Insert the Air Tube Plug into the Cuff Jack. Make 3-5 deep inhalations and exhalations before taking a measurement and relax. Do not move or speak and do not tighten your arm.

2. Press **□**.

3. All symbols will appear on the display for a short time, two short sound signals will be given and the Device will inflate automatically until the cuff is fully inflated. The figures on the screen will count back. The pulse symbol **■** will start flickering.

As the blood pressure and pulse are measured during air deflation from the arm cuff keep quiet and do not move your arm and do not tighten your arm muscles.

4. When the measurement is complete the cuff will gradually deflate. The figures on the screen will count back. The pulse symbol **■** will stop flickering.

As the blood pressure and pulse are measured during air deflation from the arm cuff keep quiet and do not move your arm and do not tighten your arm muscles.

5. After reaching the level of 190 mmHg the cuff will gradually deflate. The figures on the screen will count back. The pulse symbol **■** will stop flickering.

As the blood pressure and pulse are measured during air deflation from the arm cuff keep quiet and do not move your arm and do not tighten your arm muscles.

6. After reaching the level of 190 mmHg the cuff will gradually deflate. The figures on the screen will count back. The pulse symbol **■** will stop flickering.

As the blood pressure and pulse are measured during air deflation from the arm cuff keep quiet and do not move your arm and do not tighten your arm muscles.

7. After reaching the level of 190 mmHg the cuff will gradually deflate. The figures on the screen will count back. The pulse symbol **■** will stop flickering.

As the blood pressure and pulse are measured during air deflation from the arm cuff keep quiet and do not move your arm and do not tighten your arm muscles.

8. After reaching the level of 190 mmHg the cuff will gradually deflate. The figures on the screen will count back. The pulse symbol **■** will stop flickering.

As the blood pressure and pulse are measured during air deflation from the arm cuff keep quiet and do not move your arm and do not tighten your arm muscles.

9. After reaching the level of 190 mmHg the cuff will gradually deflate. The figures on the screen will count back. The pulse symbol **■** will stop flickering.

As the blood pressure and pulse are measured during air deflation from the arm cuff keep quiet and do not move your arm and do not tighten your arm muscles.

10. After reaching the level of 190 mmHg the cuff will gradually deflate. The figures on the screen will count back. The pulse symbol **■** will stop flickering.

