

Little Doctor

Digital Blood Pressure Monitor Instruction Manual

Ciśnieniomierz elektroniczny automatyczny do pomiaru ciśnienia tętniczego krwi i pulsu Instrukcja obsługi



ENG

POL

TABLE OF CONTENTS

PARTS AND COMPONENTS
GENERAL INFORMATION
INDICATIONS FOR USE 4
OPERATION PRINCIPLE 4
APPLIED NEW TECHNOLOGIES LD
RECOMMENDATIONS ON CORRECT MEASUREMENTS
POWER SUPPLY
BATTERY INSTALLATION 6
USE OF THE DEVICE WITH THE POWER SOURCE
SETTING OF DATE AND TIME 7
CORRECT POSITION DURING MEASUREMENT7
CUFF PREPARATION
MEASUREMENT PROCEDURE
FORCED DEFLATION FROM A CUFF10
MEMORY FUNCTION 10
MEMORY CLEARING11
INFORMATION ABOUT ERRORS11
CARE, STORAGE, REPAIR AND DISPOSAL
TROUBLESHOOTING TIPS 12
WARRANTY
TECHNICAL SPECIFICATIONS 14
CERTIFICATION AND STATE REGISTRATION 15
INFORMATION IN ENGLISH 2
INFORMATION IN POLISH

PARTS AND COMPONENTS



- 1. Electronic Unit.
- 2. LCD.
- 3. Arm cuff jack.
- 4. Air tube with air plug.
- 5. Arm cuff CUFF-LDA.
- 6. Instruction Manual.
- 7. Warranty card.

- 8. Electrical power supply LD-N057.
- 9. Power source Jack.
- 10. M1 (memory 1).
- 11. M2 (memory 2).
- 12. () (Power ON/OFF).
- 13. Storage Case.
- 14. Power elements.

GENERAL INFORMATION

This Instruction Manual is designed to assist the user with safe and effective operation of the automatic digital Device for measurement of blood pressure and heartbeat rate LD, modification LD3 (LD3a) (hereinafter – the "Device"). Use this Device according to the rules described in this Manual. Operate the Device only as intended. Do not use the Device for any other purposes. Read and understand the whole Instruction Manual, in particular "Recommendations on Correct Measurement".

INDICATIONS FOR USE

Use this Device to measure your systolic and diastolic blood pressure and heartbeat rate in patients aged from 15. This Device is recommended for use by persons with unstable blood pressure or known arterial hypertension at home as an addition to medical surveillance. The cuff is designed for the upper arm with the circumference approximately from 25 to 36 cm.

OPERATION PRINCIPLE

This Device uses the oscillometric method of blood pressure and pulse rate measurement. Wrap the cuff around your upper arm and it starts to be inflated automatically. The sensitive element of the Device feels the weak pressure oscillations in the cuff generated by widening and contraction of the brachial artery in response to every heartbeat. Pumping in is ceased when cuff is adequately pumped in to determine diastolic and systolic pressure (the amplitude of the pressure waves is measured, converted into millimeters Hg and shown on the display as figures) after which air is released from cuff. The Device has 2 memories, by 30 cells, for storage of measurement results. Remember that the Device will not maintain the mentioned accuracy of a measurement if it is used or stored at a temperature or humidity other than those specified in Technical Specifications of this Manual. We are warning about possibility of mistakes in blood pressure measurement with this Device in persons with pronounced cardiac arrhythmia. Consult the doctor concerning blood pressure measurement of your child.

APPLIED NEW LD TECHNOLOGIES



Inflation Measuring System (IMS) – Algorithm enabling to measure pressure in the process of cuff pumping in.



Indication of arrhythmia – special symbol «**Q**» on Device display informs about availability of irregular pulse; in this case.



Scale WHO – classification of measurement results according to recommendation of World Health Organization (WHO).



Fuzzy Algorithm – is the algorithm for processing the measurement values with regard to peculiarities of the man's heartbeat, thus, ensuring high measurement accuracy.

WARNING! This Device may be used only with cuff Cuff-LDA, size 25-36 cm (delivered in a set with the Device).

RECOMMENDATIONS ON CORRECT MEASUREMENTS

1. For correct measurement you should know that THE BLOOD PRESSURE IS SUBJECT TO SHARP VARIATIÓNS EVEN WITHIN THE SHORT TIME INTERVALS. The blood pressure depends on many factors. It is usually lower in summer and higher in winter. The blood pressure varies together with the atmospheric pressure, depends on physical loads, emotional excitement, stresses and dietary regime. Drugs, drinking alcohol and smoking produce significant effect. Even the very procedure of blood pressure measurement in a polyclinic sends the blood pressure high in many people, thus, the blood pressure measured at home often differs from the values received in a polyclinic. As the blood pressure tends to rise at low temperatures, make measurements at an indoor temperature (approximately 20° C). If this Device stayed under a low temperature, keep it for at least 1 hour at an indoor temperature before use, otherwise the measurement result may be incorrect. During a day the difference in readings for healthy people may be 30-50 mmHg of systolic pressure and to 10 mmHg of diastolic pressure. The dependence of the blood pressure on various factors is individual for each person. Accordingly, it is recommended to keep a special book with blood pressure records. ONLY A CERTIFIED DOCTOR USING YOUR RECORDS IS CAPABLE TO ANALYZE THE TENDENCY OF YOUR BLOOD PRESSURE VARIATIONS.

2. At cardiovascular and some other diseases requiring blood pressure monitoring make measurements in the hours fixed by your attending doctor. REMEMBER THAT THE DIAGNOSTIC AND ANY TREATMENT OF HYPERTENSION MAY BE CONDUCTED ONLY BY A CERTIFIED DOCTOR ON THE BASIS OF BLOOD PRESSURE VALUES OBTAINED BY THIS DOCTOR. TAKING OF DRUGS AND THEIR DOSES SHOULD BE PRESCRIBED ONLY BY YOUR ATTENDING DOCTOR.



3. At such disorders as deep vascular sclerosis, weak pulse wave and also in patients with the prominent distortions of cardiac rhythm it may be difficult to measure the blood pressure accurately. IN SUCH CASES CONSULT A CERTIFIED DOCTOR ABOUT APPLICATION OF THE ELECTRONIC DEVICE.

4. KEEP QUIET DURING A MEASUREMENT TO OBTAIN THE ACCURATE VALUES OF YOUR BLOOD PRESSURE WITH THE ELECTRONIC DEVICE. Measure your blood pressure in the calm and comfortable conditions at the indoor temperature. No eating an hour before measurement; no smoking, taking tonic agents, alcohol 1.5-2 hours before measurement.

5. The accuracy of blood pressure measurement depends on whether the cuff matches the size of your arm. THE CUFF SHOULD NOT BE TOO SMALL OR TOO LARGE.

6. Wait 3 minutes between measurements for the blood to restore its circulation. However, the persons with prominent atherosclerosis due to considerable loss of vascular elasticity may need to increase the wait time between measurements (10-15 min.). This also refers to the patients suffering for long from diabetes.

For more accurate determination of blood pressure it is recommended to make a series of 3 consecutive measurements and to use the average value.

POWER SUPPLY OF THE DEVICE

BATTERY INSTALLATION

1. 1. Open the cover of the battery compartment and install 4 "AA" size batteries according to polarity marked inside the compartment. Do not use much force to remove the cover of the battery compartment (fig. 1).

- 2. Close the battery cover.
- Replace all batteries when the Low Battery Indicator " ", appears on the screen or when there is no any indication on the screen. The Low Battery Indicator does not show the discharge level.



Fig. 2

- The batteries supplied with the Device are intended for check of the Device performance at sale and their service life may be shorter than of the recommended batteries.
- Replace all four batteries at the same time. Do not use the waste batteries.
- If the Device is unused for a long time, remove all batteries.
- Do not leave the waste batteries in the Device.
- Rechargeable power cells, type AA, may also be used.

USE OF THE DEVICE WITH THE POWER SOURCE

The manufacturer recommends application of the stabilized power source LD-N057 (it is attached to modification LD3a).

The jack for connection of the power source is on the right side of the Device.

To use the device with electric power source (EPS), connect plug connector of EPS to the Device and insert plug of EPS into mains socket.

Measurement over, switch off the Device having depressed button, remove plug of EPS from mains socket and disconnect plug connector from the Device.

SETTING OF DATE AND TIME

To shift to mode of date and time setting, it is necessary, holding M1, depress (). Selected parameter will blink. Change of selected parameter toward increasing is taking place when depressing M2, towards decreasing – when depressing M1. To shift to setting of the next parameter, year/month/hours/min, it is necessary to depress ().

If no actions are done in mode of date and time setting for more than 1 minute, the Device is independently changed over to mode of date and time indication.

When power cells are replaced, date and time will zero.

Measurement of pressure and pulse is possible without date and time setting. Indication of date and time is not switched off. Current consumption in this case is negligible.

CORRECT POSITION DURING MEASUREMENT

- 1. Sit at a table so that during blood pressure measurement your hand rests on its surface. Be sure that the cuff is placed approximately at the level of your heart and that your arm lies freely on the table and does not move (fig. 3).
- 2. You can measure the blood pressure lying on the back. Look at the ceiling, keep quiet and do not move during measurement. Be sure that the cuff is placed approximately at the level of your heart (fig. 4).







CUFF PREPARATION

1. Insert the cuff end for about 5 cm into a metal ring as shown in the figure 5.

2. Apply the cuff to your left upper arm so that the air tube is directed to your palm. If the measurement on your left arm is difficult, you may use your right arm. In this case remember that the readings may differ by 5-10 mmHg and even more (fig. 6).

3. Wrap the cuff around your upper arm so that the bottom of the cuff is approximately 2-3 cm above your elbow. The sign "ARTERY" should be over the arm artery (fig. 7).

4. Fix the cuff so that it fits tightly to the arm, but see that it is not overtight. Too tight or too free placement of the cuff may give inaccurate readings (fig. 8).

5. On the fixed cuff the sign «INDEX» should point to the area «NORMAL (25-36 cm)». It means that the cuff is chosen correctly and fits the size of your upper arm. if the sign points to the area marked «**411**» the cuff is too small and the readings will be higher. If the sign points to the area marked **411**» the cuff is too large and the readings will be lower (fig. 9).



Fig. 5





Fig. 7





6. If the arm has a conic form, the cuff should be put on with a spiral movement as shown in the figure (fig. 10).

7. If the rolled-up sleeve squeezes the arm interfering with free blood flow the Device may give inaccurate figures not corresponding to your actual blood pressure (fig. 11).

MEASUREMENT PROCEDURE

1.Insert the Air Tube Plug into the Cuff Jack. Make 3-5 deep inhales and exhales before taking a measurement and relax. Do not move, do not speak and do not toughen your arm.

2. Press 🕛.

3. All symbols will appear on the display screen for a short time, two short sound signals will be given and the Device will inflate automatically the air into the cuff (fig. 12),

4. After level 18-20 mm Hg is reached, cuff will start to be pumped in, first slowly, then quickly. During pumping in sign «♥» is blinking. Algorithm IMS makes it possible to determine systolic and diastolic pressure in the process of pumping in (fig. 13).

5. At the end of measurement, depending on classification of result according to scale of World Health Organization (WHO), short sound signals will ring out: 1 signal – normal pressure; 2 signals – increased pressure within the norm; 3 signals – hypertension. After that the device will release all air from cuff and measured arterial blood pressure, pulse rate and time of measurement will be displayed on the display (fig. 14). Symbol «fig.», is blinking on display, reminding that to retain results,

it is necessary to choose memory 1 or 2, having depressed M1 or M2, respectively. If not to choose memory within 3 minutes, the result will not be memorized, and the device is automatically switched over to mode of date and time indication.



Fig. 10











ENG

If irregular pulse rhythm is detected during measurement, symbol of arrhythmia «♥», will appear and blink upon measurement end. Appearance of arrhythmia indicator may also be caused by movement of body during measurement. During periodical appearance of this indication «♥» apply to You attending doctor.

Apart from numerical value of pressure, result is also displayed on scale WHO (fig. 14). Scale WHO – three-color scale of classification of received value of arterial pressure, according to recommendation of World Health Organization (fig. 15). The scale is located in the left corner of display.

6. Press () to shift to mode of date and time setting.

TO OBTAIN THE ACCURATE RESULT MAKE INTERVAL BETWEEN MEASUREMENTS TO RESTORE THE BLOOD CIRCULATION. WAIT FOR AT LEAST 3 MINUTES BEFORE MAKING A NEW MEASUREMENT.

THE DATA WILL BE KEPT IN THE MEMORY EVEN IF THE DEVICE IS STORED WITHOUT BATTERIES. TO DELETE ALL VALUES STORED IN THE MEMORY YOU SHOULD MAKE ACTIONS DESCRIBED IN "MEMORY FUNCTION".

If no actions are done in mode of date and time setting for more than 3 minute, the Device is independently changed over to mode of date and time indication.

FORCED DEFLATION FROM A CUFF

For rapid air release from of the arm cuff during arm cuff inflation or during a measurement (slow deflation) press the \bigcirc Button. The arm cuff quickly completely deflates and the Device is independently changed over to mode of date and time indication.

MEMORY FUNCTION

1. Result of each measurement (pressure, pulse, time and date of measurement) may be kept in the Device memory. For this purpose, after measurement, within not more than 3 minutes, memory M1 or M2 shall be selected for memory storage.

IF THE NOTICE ON ERROR APPEARS THE MEASUREMENT RESULT WILL NOT BE STORED.

2. Up to 30 measurement results and mean value of last 3 measurements may be kept in each Device memory. When the number of measurements exceeds 30, the oldest record is deleted to save the most recent values.

3. Press the M1 (or M2) to see the figures stored in the memory. At the first depression of M1 (or M2) mean value of 3 last measurements, kept in memory M1 (or M2), will appear on the screen, Mean value of pulse rate is displayed by symbol "A" for 3 last measurements. At repeated depression of M1 (or M2) indicator of selected memory M1 (or M2) «



 $(\operatorname{ore}_{A}^{T_{p}})$ and number of memory cell will appear on the screen, and in 1 second its content is displayed (fig. 17). When content of memory cell is displayed, date and time of measurement are displayed alternately in the display lower line. Each depression of M1 (or M2) causes shifting to the next memory cell. If there are no saved results in memory, at first depression of M1 (or M2) current date and time are displayed on display Device.

MEMORY CLEARING

To delete all values stored in the memory, press the M and hold it down for more than 5 seconds. Symbols "CIr" will appear on the screen and all values will be deleted from the memory (fig. 18).



ENG



INFORMATION ABOUT ERRORS

Indication	Likely cause	Methods of correction
Err	The arm cuff is applied incorrectly or the air tube plus is inserted not tightly enough.	Be sure that the arm cuff is applied correctly and the plug is inserted tightly. Repeat the whole measurement procedure.
	Measurements cannot be made due to hand movement or talking during measurements.	Repeat the measurement fol- lowing strictly the recommen- dations of this Manual.
	Batteries are discharged.	Replace the batteries for new ones.

CARE, STORAGE, REPAIR AND DISPOSAL

- 1. Keep this Device from exposure to higher humidity, direct sunlight, shocks, vibration. THIS DEVICE IS NOT WATERTIGHT.
- 2. Do not keep and use this Device near heating installations and open fire.
- 3. If the Device was stored at a temperature below the freezing point, keep it at least for 1 hour in some warm place before use.
- 4. Remove the batteries if the Device will be unused for a long time. Battery leaking may damage the Device. KEEP BATTERIES OUT OF REACH OF CHILDREN!
- 5. Keep the Device clean and protect it from dust. Use the dry soft cloth to clean the Device.

- 6. Keep the Device and its components away from water, solvents, spirit and petrol.
- Protect the arm cuff from contacting on sharp things; do not stretch or fold tightly the arm cuff.
- 8. Do not subject the Device to strong shocks, such as dropping on the floor.
- 9. This Device does not contain special controls to adjust the measurement accuracy. It is prohibited to open individually the electronic block. Repair the Device only in authorized organizations.
- 10. On expiration of the warranted service life apply from time to time to authorized repair organizations to check the technical condition of the Device.
- 11. Dispose of the Device and its components according to the application local regulations. No special requirements to disposal of this Device are defined by the manufacturer.
- 12. The arm cuff may withstand multiple sanitary treatments. The internal tissue surface of the arm cuff (contacting on arm) may be cleaned with cotton ball soaked in 3%-solution of hydrogen peroxide. At long use the partial color fading of the tissue coating of the arm cuff is possible. Washing and ironing of the arm cuff are not allowed.
- 13. Before use of EPS, check integrity of power supply cord of EPS.

L.	ROODLESHOOTING	111.5
PROBLEM	LIKELY CAUSE	METHOD OF CORRECTION
After pressing the O/I no indication on the display.	Discharge of batteries. Polarity of batteries is not observed. Battery terminals are contaminated. The power supply is not inserted into a power outlet.	Replace all batteries for new ones. Install batteries correctly. Clean the terminals with dry cloth. Insert the power supply into the power outlet.
Every time the blood pressure is different.	Check that the arm cuff is level with your heart?	Take the correct position for measurement.
low/high.	Perhaps you talk or move your arm during the measurement?	Calm down before the measurement. Keep silence and quiet during measurement.

TROUBLESHOOTING TIPS

ENG

PROBLEM	LIKELY CAUSE	METHOD OF CORRECTION
Measurements of the pulse rate are too high/low.	Perhaps you talk or move your arm during the measurement? Perhaps you make measurement directly after physical load?	Keep silence and quiet during measurement. Repeat the measurement at least in 5 min.
Impossibility to make a large of number of measurements.	Application of poor batteries.	Use only alkali batteries of well-known manufacturers.
Independent changing over to mode of date and time.	System of automatic change-over to mode of date and time indication operates.	This is not a defect. This is not a trouble. The Device is automatically changing over to mode of date and time indication in 3 minutes after use.

ENG

If regardless of the above recommendation you are unable to get correct measurement results, stop the use of this Device and apply to a maintenance organization (addresses and telephones of authorized organizations may be found in the warranty card). Do not try to adjust the internal mechanism by yourself.

WARRANTY

- 1. The following LD product is covered by warranty for the period specified in the warranty card.
- 2. The warranty liabilities are contained in the warranty card given at the sale of this Device to a purchaser.
- 3. The addresses of organizations for warranty maintenance are given in the warranty card.

TECHNICAL SPECIFICATIONS

Measurement method	oscillometric with Fuzzy Algorithm and Inflation Measuring System
Display	LCD, four-line display
Pressure indication range in an arm cuff, mmHg	from 0 to 300
Measurement range: pressure in an arm cuff, mmHg pulse rate, 1/min	from 40 to 260 from 40 to 160
Range of admissible error: of air pressure in an arm cuff, mmHg at pulse rate measurement, %	±3 ±5
Inflation	automatic (air pump)
Deflation	automatic
Date and time	yes
Memory	2x30 recent measurements + average value of the last three measurements
Power source, V	6
Type of power supply	4 "AA" size batteries (LR6) or adapter not less than 600 mA
Max power intake, W	3,6
ADAPTER LD-N057	
Output voltage,V Max load current, mA Input voltage Dimensions , mm Weight, kg Plug: Polarity of terminals Internal diameter, mm External diameter, mm Length of plug contact, mm	$6 \pm 5\%$ not less than 600 ~200-240 V, 50/60 Hz $64 \times 70 \times 43$ not more 0,3 «-» internal 2.1 ± 0.1 5.5 ± 0.1 10 ± 0.5
Operation conditions: Temperature, ℃ Relative humidity, % Rh	from 10 to 40 85 and lower
Storage and transportation conditions: Temperature, °C Relative humidity, % Rh	from -20 to 50 85 and lower

Cuff size:	larger for adults (upper arm circumference 25-36 cm)
Dimensions: Size (electronic block), mm Weight (without package, case, batteries and adapter), g	131 x 146 x 130 424
Completeness	electronic block, Cuff-LDA (in a set with a tube and plug), 4 batteries, power source LD-N057, instruction Manual, warranty card, a case, package
Year of manufacture	year the manufacture is given in the bottom of the Device body in a serial number after symbols "AA"
Symbol definition	★ Type BF applied part. I Read Instruction Manual. C€0123 European Union Approval.

ENG

Technical characteristics may be changed without preliminary notification to improve the operation and quality of the product.

CERTIFICATION AND STATE REGISTRATION

This unit manufacturing is certified according to international standard ISO 13485:2003. Units LD6 comply with the requirements of European Directive MDD 93/42/EEC, international standards, EN980, EN1041, EN1060-1, EN1060-3, EN10601-1-2, ISO 14971, EMC (IEC 60601-1-2:2001/A1:2004, CISPR 11:2003/A2:2006, GROUP 1, CLASS B, IEC 61000-3-2:2005, IEC 61000-3-3:1994/A2:2005.

Power source LD-N057 complies with international standard EN 55022 Class A, protection level against electric shock: Class II, Type B.

☑ Complaints and requests should be addressed to:

Little Doctor Europe Sp. z o.o. 57G Zawila Street, 30-390, Krakow, Poland Service phone: +48 12 2684748, 2684749

Manufactured under control:

Little Doctor International (S) Pte. Ltd., 35 Selegie Road # 09-02 Parkline Shopping Centre, Singapore 188307

Manufacturer:

Little Doctor Electronic (Nantong) Co. Ltd., No.8, Tongxing Road Economic & Technical Development Area, Nantong 226010, Jiangsu, PEOPLE'S REPUBLIC OF CHINA

Distributor in Europe:

Little Doctor Europe Sp. z o.o. 57G Zawila Street, 30-390, Kraków, Poland Sales Office phone: +48 12 2684746, 12 2684747, fax: +48 12 268 47 53 E-mail: biuro@littledoctor.pl www.LittleDoctor.pl



 Omega
 Authorized Representative in the EU:

 Shanghai International Trading Corp.GmbH
 Eiffestrasse 80, 20537 Hamburg, Germany.

www.LittleDoctor.pl

EC



LITTLE DOCTOR INTERNATIONAL (S) PTE. LTD. Yishun Central P.O. Box 9293 Singapore 917699, Fax: 65-62342197, E-mail: Id@singaporemail.com

> REP Shanghai International Trading Corp.GmbH Eiffestrasse 80, 20537 Hamburg, Germany. Tel: 0049-40-2513175

Registered trade marks of Little Doctor International (S) Pte. Ltd.
Little Doctor International (S) Pte. Ltd., 2007-2015

E347/1504/2